

Hygienic Lanyards for Audio Tour Equipment

Silicone Lanyards – Opus+

Easily disinfected with 70% Isopropyl Alcohol or 3% Hydrogen Peroxide and 97% Water Solution

Multiple color options

(not available for Opus Click or Touch)



Plastic Lanyards – Opus Click & Opus Touch

Easily disinfected with 70% Isopropyl Alcohol or 3% Hydrogen Peroxide and 97% Water Solution

Multiple color options

(not available for Opus+)



Plastic Lanyards – Opus+ & SmartScreen Devices

Easily disinfected with 70% Isopropyl Alcohol or 3% Hydrogen Peroxide and 97% Water Solution

Multiple color options and can be custom branded for additional fee

(not available for Opus Click or Touch)



TPU Plastic Sleeves for Audio Tour Equipment




Plastic Equipment Sleeves

TPU (thermoplastic polyurethane) Plastic –can be disinfected between uses with a 70% Isopropyl Alcohol or 3% Hydrogen Peroxide cleaning solution
Magnetic closure. Can be used with speaker or earbuds.

Features:

- Antibacterial coating
- Optional Built-in External Earphone Jack (2.5mm or 3.5mm)



Hygienic Accessories for Audio Tour Equipment

Single Use Earbuds

Individually packaged earbuds
Available with 2.5mm and
3.5mm jacks



Plastic Equipment Bag

Single use, disposable, plastic

Please contact for pricing

acoustiguide

sales@acoustiguide.com
www.acoustiguide.com